



Disposable Core Biopsy System SABD™




Helping Answer Life's
Most Important Questions®




| Features |
|-----------------------------------|
| Full specimen notch exposure |
| Bevel or Trocar tip configuration |
| Echogenic tip |
| Etched depth marks |
| ISO color coded |
| Light weight ergonomic design |

Traditional Bevel Tips
are available in every gauge
and length combination.



Trocar Tips provide
greater cutting power
for dense tissue.



| Benefits |
|---|
| Provides full sampling, easy specimen access |
| Provide a choice for the user when sampling in dense tissue |
| Enhances visibility under ultrasound |
| Provides confirmation of position |
| Provides easy identification |
| Provides comfort, ease of use for CT applications |

| SABD Product No. | Description | Per Box |
|------------------|--------------------------------|---------|
| SABD-1409-15 | SABD 14gx 09cm,20mmT bevel tip | 10 |
| SABD-1409-15-T | SABD 14gx 09cm,20mmT bevel tip | 10 |
| SABD-1415-15-T | SABD 14gx15cm,20mmT bevel tip | 10 |
| SABD-1609-15 | SABD 16gx 09cm,20mmT bevel tip | 10 |
| SABD-1609-15-T | SABD 16gx 09cm,20mmT bevel tip | 10 |
| SABD-1615-15 | SABD 16gx 15cm,20mmT bevel tip | 10 |
| SABD-1809-15 | SABD 18gx 09cm,20mmT bevel tip | 10 |
| SABD-1809-15-T | SABD 18gx 09cm,20mmT bevel tip | 10 |
| SABD-1815-15 | SABD 18gx 15cm,20mmT bevel tip | 10 |
| SABD-1815-15-T | SABD 18gx 15cm,20mmT bevel tip | 10 |
| SABD-1820-15 | SABD 18gx20cm,20mmT bevel tip | 10 |
| SABD-1820-15-T | SABD 18gx 20cm,20mmT bevel tip | 10 |
| SABD-2006-10-T | SABD 20gx 06cm,15mmT bevel tip | 10 |
| SABD-2009-15-T | SABD 20gx 09cm,20mmT bevel tip | 10 |
| SABD-2015-15 | SABD 20gx 15cm,20mmT bevel tip | 10 |
| SABD-2020-15-T | SABD 20gx 20cm,20mmT bevel tip | 10 |

The Products, photographs, images and words belong to their respective companies.